



Suite 3.21 Innovation Central

10 John Williams Boulevard South, Central Park, Darlington, DL1 1BF
+44 (0) 1325 353776

support@pukka-j.com www.pukka-j.com

Release Notes

Pukka-J DICOM Explorer 6.4.6

(Released 10/10/2023)

Patient Merge Issue

Fixed an issue with merging patients in the RIS.

Deletion Speed

Improved deletion speed of large studies.

RIS HL7 Reports

Modified the ORU handling in the RIS by replacing the \.br\ characters with carriage returns.