



Pukka-j, Blackwell Farm  
Bishops Cannings  
Devizes, Wiltshire  
SN10 2JZ, UK

Tel: +44 (0) 1380 860044  
Fax: +44 (0) 1380 860066  
Email: [support@pukka-j.com](mailto:support@pukka-j.com)  
Web: [www.pukka-j.com](http://www.pukka-j.com)

## Release Notes

### **mint.jar 1.22.07**

(Released 06/01/2021)

#### **JPEG2000 Node**

A node dedicated to the handling of JPEG2000 data has been added. This provides the option of using an external decompression executable process to do the decompression by writing a temporary cache file and a library option to perform the decompression using a CODEC library.

The Kakadu implementation of the JPEG2000 codec has been added as a library option.

#### **DicomFileOutput Node**

Fixed an issue with a UID check on currently active files being written by the output node.

#### **DicomOutputPort Node**

In the event that an illegal file is received due to the pixel data being an odd length, the output port will fix the issue by increasing the tag length by one and adding a padding byte to the last frame in the image data tag.

#### **Stream Splitter Node**

Fixed a splitter issue with the new protocol filter in that the protocol filter was masking the AET and tag checking.

#### **Temp Cache / Forwarder Node**

Added a file exist check to the forward job because if the same study was sent, it could create multiple records to forward but it will only be done once then deleted which created a file reference error and left the forwarder in a loop trying to forward a file that was already deleted in a previous forward job.