



Pukka-j, Blackwell Farm  
Bishops Cannings  
Devizes, Wiltshire  
SN10 2JZ, UK

Tel: +44 (0) 1380 860044  
Fax: +44 (0) 1380 860066  
Email: [support@pukka-j.com](mailto:support@pukka-j.com)  
Web: [www.pukka-j.com](http://www.pukka-j.com)

## Release Notes

**cloe.js.gz**

**nexus.jar 2.18.02**

**mint.jar 1.20.03**

**pacsoology.jar 1.16.02**

**pjk.jar 1.15.02**

(Released 27/01/2020)

### **DICOM to HL7 Node**

A connection closed issue with the HL7 output port was resolved.

### **Implemented DICOM Patient Level Move SOP Class**

The DICOM Input Port was extended to allow the Patient Root Move protocol. This option is now available in the user interface of the Dicom Stream Input Node.

See individual library release notes for more details.