



Pukka-j, Blackwell Farm
Bishops Cannings
Devizes, Wiltshire
SN10 2JZ, UK

Tel: +44 (0) 1380 860044
Fax: +44 (0) 1380 860066
Email: support@pukka-j.com
Web: www.pukka-j.com

Release Notes

pjk.jar 1.15.10

(Released 05/08/2020)

PID Issuer

Added some control over how the issuer on the patient ID is handled in the DICOM Query Result.

DICOM Dose SR

Extended the DicomBasicTextSR to create a DicomDoseSR for holding PJ Dose and post processing information.

Time Stamp Authority

Updated the timestamp authority used when signing the library.

tsaurl=<http://sha256timestamp.ws.symantec.com/sha256/timestamp>

Encapsulated PDF Series Date and Time

Added series date and time when creating an encapsulated PDF object.

Logging

Updated some logging output.