



Pukka-j, Blackwell Farm
Bishops Cannings
Devizes, Wiltshire
SN10 2JZ, UK

Tel: +44 (0) 1380 860044
Fax: +44 (0) 1380 860066
Email: support@pukka-j.com
Web: www.pukka-j.com

Release Notes

pjk.jar 1.15.11

(Released 21/10/20)

DICOM Stream Closing

When reading a DICOM file via an `InputStream`, the DICOM class to which the stream was passed to automatically called the close method on the stream when finished. Under normal situations, this is expected, but in order to extend the classes to handle a set of DICOM files in a zip archive, where the `ZipInputStream` is the same `InputStream` for all `ZipEntries`, the `InputStream` shouldn't be closed, so an option was added when reading a DICOM File as a stream to prevent the automatic closing of the `InputStream`.

ZIP Testing

Added the byte signatures of some common ZIP files to the `OIUtil` class to test the first few bytes of a file to establish if it is a ZIP archive (or not).

JSON Formatting

There was an erroneous character added to the JSON output when adding escapes to a carriage return when formatting an output for a JSON HTTP response.